
SOLTEC

LabNews

SOLTEC and l'IRCCS Mario Negri Institute for Pharmacological Research in Milan

A partnership for innovation
and participatory research

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www.soltec.eu
e-shop

Follow us: [f /soltec.ultrasonics](https://www.facebook.com/soltec.ultrasonics)

Follow us: [in /company/soltec-srl](https://www.linkedin.com/company/soltec-srl)

Stefano Guidi President SOLTEC srl
Stefano Guidi Presidente SOLTEC srl

“Help research to help life!”

This is the slogan which underlies the tasks carried out by the research team at the Mario Negri Institute of Pharmacological Research in Milan. This private nonprofit foundation is run by Prof. Silvio Garattini, and its mission is to serve the sick and defend human health and life. In order to maintain the quality of its scientific research, the Institute has maintained its independence, siding with the sick; the Institute chooses not to patent the results of its research, but makes them freely available to the scientific community and patients. As a concrete show of support for this philosophy, in early 2015, here at SOLTEC we decided to donate some of our SONICA ultrasound equipment [ultrasonic baths and the special SONICA Ultrasonic Extractor] to the Mario Negri Institute to enable their research teams at the Environmental Chemistry and Toxicology laboratories to carry out sonication processes. This helped the researchers to extract pollutants from environmental matrices and PM10-type particulate matter from membranes. This technical and scientific magazine and the video interviews with Prof. Silvio Garattini and PhD researcher Diego Baderna on our corporate website www.soltec.it describe the ultrasonic extraction methods used at the Mario Negri Institute.

But SOLTEC's work and commitment to research and innovation does not end there!

In the early months of 2015, we developed and engineered a new system for lifting and moving washing loads, aimed at small and medium-sized enterprises in the industrial sector. This new system, which will be marketed by

SOLTEC later this year and distributed via our dealers in Italy and elsewhere, grew out of our designers' burning desire to make manual work more efficient and functional. Indeed, this system is based on a study of musculoskeletal disorders and dysfunctions, often associated with working conditions involving manual lifting and moving of heavy loads.

The innovative SOLTEC system, which is being presented in the SOLTEC LabNews journal, provides a so-called “assisted” lifting mechanism which can eliminate any strain on the part of the worker.

This system can also be combined with a practical motorized oscillator for the basket inside the ultrasonic bath, thus reducing cleaning times and greatly improving quality.

In this special issue dedicated to research and innovation, as well as the above editorial contributions, we are delighted to include other special articles, by Dr. Christian Durante from the Laboratory of Chemical Sciences at Padova University, by researchers from TUM (Technische Universität München) and many others. The magazine would not have been complete without a brief presentation of our recent trade fair experiences in places as far afield as MEDICA in Düsseldorf, ADF in Paris, IDS in Cologne, ArabLab in Dubai andACHEMA in Frankfurt, which have put us in touch with new customers in Europe, Asia and Arab countries.

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ANALYSIS OF A NEW LOAD LIFTING AND HANDLING SYSTEM MANUFACTURED BY SOLTEC

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Musculoskeletal diseases and disorders are becoming increasingly important in the fields of occupational medicine and prevention.

Over the years, a clear association has been observed between such disorders and working conditions prevailing in various companies, mostly related to tasks defined by or connected with certain manual operations, which can be physically very demanding. For that reason, such disorders are widespread, which has repercussions in terms of economic and social costs.

Moreover, 2012 data from INAIL show that work-related musculoskeletal disorders are on the increase, accounting for about 56% of all occupational diseases recognized in Europe; in Italy, in 2012, there were 26,000 reported cases of work-related diseases related to bone and joint problems.

This fact in itself highlights just how much we tend to underestimate the predisposing factors for these diseases when assessing individual work stations and, as a result, during the risk assessment process as a whole.

According to analyses by the European Agency for Safety and Health at Work, these disorders seem to be caused mainly by manual handling of loads, by heavy physical work, by unsuitable working postures, by

1

Figure 1 - New modular intensive cleaning system for industry service and maintenance

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Figure 2 - New modular system with lifting device and automatic oscillating movement of the basket

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2

performing repeated, sometimes strenuous, movements of the upper limbs and, as an aggravating factor, also by the vibrations produced by work equipment or vehicles.

Annex XXXIII to Italian Legislative Decree No 81/2008 identifies the risk factors that the employer has to watch out for and which may be involved in manual handling of loads, during specific work phases.

These tasks definitely need to be assessed so that preventive measures can be taken to mitigate them, or to eliminate them from individual work phases.

Soltec designs and manufactures ultrasonic cleaners for professional use in the chemical, pharmaceutical, medical, electronics, military and aerospace industries. Given the wide variety of uses that such cleaners can be put to, there are different sizes and workloads, which means that the size and weight of the article to be cleaned will also change.

Analyzing the work flow associated with the machine, it appears evident that, especially when loading and unloading workpieces, where the loads exceed 5-10 kg or more, the operator performs lifting and handling movements that

in particular strain the upper limbs and spine, especially when carried out on a regular basis throughout the working day.

On the basis of these mandatory indications and with a view to making manual labour more efficient and functional, Soltec decided to engineer a lifting mechanism for lifting baskets; this mechanism, with an adjustable counterweight, makes it possible to perform what we might call "assisted" lifting, thus limiting the effort made by the operator or even totally removing it.

From an ergonomic and research standpoint, we initially analyzed the movements of the user, verifying what the problematic movements and operations were in terms of work and lifting ergonomics.

Thereafter, on the basis of measurements taken at the workstation and from applicable technical standards* data using, amongst other techniques, high-speed cameras and digital crane scales, we were able to analyze the efforts made by the operator for given load ranges and consequently develop a product that was both suitable for making the workstation ergonomically favourable and functional in terms of productivity and load management.


3

Figure 3 - New modular system with automatic oscillating movement of the basket

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The professional lift designed by SOLTEC, with its modularity and its adjustability to the physique of the operator (height, limb size, weight to be lifted, etc.), enables the operator to lift loads of up to 30 kg keeping the spine in an upright position and the shoulders close to the body, so as not to place any strain on the muscles and thus protect him/her from the effects of muscle tension in the neck and shoulders.

In the case of Soltec's ultrasonic cleaners, at our assisted workstations the operator will place the object to be cleaned inside the basket, on a shelf at a height that can easily be adjusted in order to minimize strain, and will then raise the loaded basket using the lifting "handlebar" to the right height to place it inside the ultrasonic cleaner. After cleaning, the basket will be lifted out, in assisted mode, and placed back on the work surface, thus making it easy to remove the object.

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Franco Beltrame

Industrial Hygienist and
Occupational Safety Technician
AL.MA. Consulenza srl

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SOLTEC and IRCCS Mario Negri Institute for Pharmacological Research in Milan A partnership for innovation and participatory research

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At the Mario Negri Institute in Milan, researchers are using increasingly modern equipment and technologies.

Here at SOLTEC, we decided to provide them with our sonication equipment free of charge to help them carry out their research.

Diego Baderna, PhD student and researcher at Laboratory of Environmental Chemistry and Toxicology, uses SONICA Ultrasonic baths to extract pollutants from environmental matrices, while the special SONICA Ultrasonic Extractor is used in two European projects to remove PM10-type particulate matter from membranes.

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In the laboratory of Environmental Chemistry and Toxicology (Department of Environment Health Sciences) of "IRCCS Istituto di Ricerche Farmacologiche Mario Negri" in Milan, researchers are studying the impact of environmental pollutants on human health and on the environmental receptors. The SONICA 3200M ultrasonic bath, currently supplied to the laboratory, is a very useful tool in the preparation of soil leachates according to the ANPA (2000) and APAT (2002 and 2004) guidelines which will be then analyzed by in vitro toxicological assays on human cells or ecotoxicological assays.

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Figure 2 - Preparation of soil leachates according to the ANPA and APAT.

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The large tank of the instrument allows the simultaneous housing of more reaction chambers, generally consisting of 100 or 250 ml glass flasks, depending on the amount of sample to be extracted. In this way, in addition to the simultaneous processing of multiple samples, it is possible to obtain a comparable extraction efficiency among the samples processed in the same series, ensuring the possibility of comparison. The large tank also allows the insertion of blocks of synthetic ice to cool the temperature of the water in the bath. Using the same technique, the SONICA 3200M was also used to obtain aqueous leachates of earth and rock from mechanized tunneling, comminuted solid waste and ashes while, coupled with organic solvents (ethanol, methanol, acetone: hexane mixture), it was used for extraction of organic pollutants from waste shredding, earth and rock from excavation, especially when the use of high temperature extraction (soxhlet) is not suitable to preserve the integrity of the sample or of the analyzed pollutant (for example PAHs are sensible to temperature > 40°C).

Milan, June 30th 2014

Diego Baderna

*PhD student and junior researcher
at Laboratory of Environmental
Chemistry and Toxicology*

4

3

Figure 3 - Preparation of soil leachates according to the ANPA and APAT.

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Figure 4 - Glass flask during ultrasonic extraction with Sonica Bath.

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The SONICA 3200M ultrasonic bath, currently supplied to the laboratory, is a very useful tool in the preparation of soil leachates

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Figure 5 - Glass flask after ultrasonic extraction with Sonic a Bath

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To do this, Istituto Mario Negri uses:

SONICA® 3200 M S3 and Sonica ultrasonic extractor

TECHNICAL CHARACTERISTICS OF SONICA 3200:

Tank dimension: 300X240X100	US Peak Power: 360 W
Capacity Litres: 6	Power Heating: -
External dimension: 325X270X260 mm	Drain cock: Ø 7mm
Frequency: 40 Khz	Sweep System

Cleaning procedures in Ultrasonic Bath SONICA 3200 EP of glassware and electrochemical cell components of Electrochemical Scanning Tunneling Microscope

The electrochemical scanning tunneling microscope (EC-STM) is a microscopy technique coupled with an electrochemical technique such as cyclic voltammetry. The EC-STM employed in the laboratory of Chemical Sciences of University of Padua permits the imaging of smooth surface with atomic resolution at the solid-liquid interface; practically this technique permits to see the surface atoms in the presence of an electrolytic solution. Therefore, EC-STM allows an in-situ evaluation of the electrocatalysis processes and to follow the electrodes surface changes upon a chemical or an electrochemical reaction at atomic resolution. Dr. Christian Durante explains us how to use the SONICA Ultrasonic cleaner with Sweep System Technology in order to eliminate any contamination from the objects surface.

Tomasz Kosmala and Christian Durante Electrocatalysis and Applied Electrochemistry Group, Department of Chemical Sciences, University of Padua, Via Marzolo 1, 35131 Padova, Italy

The electrochemical scanning tunneling microscope (EC-STM) is a microscopy technique coupled with an electrochemical technique such as cyclic voltammetry. The EC-STM employed in our laboratory permits the imaging of smooth surface with atomic resolution at the solid-liquid interface; practically this technique permits to see the surface atoms in the presence of an electrolytic solution. Therefore, EC-STM allows an in-situ evaluation of the electrocatalysis processes and to follow the electrodes surface changes upon a chemical or an electrochemical reaction at atomic resolution.

In order to obtain high quality of in-situ STM images it is necessary to eliminate possible contaminations coming from the solution, electrochemical cell, glassware, tweezers, glue used for coating of the tips etc. In Figure 1 a) and b) it is reported an example of a contaminated system; the surface has been not sufficiently cleaned so that atomic resolution is obtained but contaminants are clearly visible as light spots over the iodine-modified Au(111) surface.

“ Ultrasonic cleaning uses Cavitation micro-to-nano-sized bubbles induced by high-frequency pressure (sound) waves to agitate a liquid. ”

Figure 1 - Particolare building University of Padua

The presence of contaminants precludes in most cases the imaging of the surface and to reach atomic resolution. Furthermore the presence of nanocontaminants on the scanner tips or sample can produce fake images and therefore misleading results.

On the contrary, when a deep cleaning is assured, it is possible to obtain pictures such as those reported in in Figure 1 c) and d). In the example reported in Figure 1 c) and d) (I adsorbed on Au(111) single crystal) it is possible to see the surface atoms and there are no traces of contaminations. The easiest and most efficient method to clean all crucial equipment is to use ultrasonic bath.

The equipment employed for the acquisition of the Figure c) and d) have been cleaned in Ultrasonic Bath SONICA 3200 EP.

Ultrasonic cleaning uses Cavitation micro-to-nano-sized bubbles induced by high-frequency pressure (sound) waves to agitate a liquid. Stimulation induces a high force on the contaminants adhering to substrate such as metals, plastics, glass, rubber and ceramics. Due to the inherently small size of the bubbles and their relatively large energy ultrasonic cleaner is able to reach into small crevices and remove entrapped soils very effectively.


Figure 2 - In situ STM images obtained from a Au(111) surface in 5mM H₂SO₄ + 1mM KI solution; image a) and b) represent an example of contaminated system, c) and d) high-resolution STM image of I/Au(111) surface without trace of contamination; images parameters: a) 28nm x 28nm, I_t = 1nA, U_b = -3mV, E = -500mV vs Pt/PtI; b) 11nm x 11nm, I_t = 1nA, U_b = -3mV, E = -500mV vs Pt/PtI; c) 11nm x 11nm, I_t = 1nA, U_b = -15mV, E = -400mV vs Pt/PtI; d) 11nm x 11nm, I_t = 1nA, U_b = -15mV, E = -400mV vs Pt/PtI.

Figure 3 - Cleaning in Ultrasonic Bath SONICA 3200 EP: a) Ultrasonic Bath SONICA 3200 EP, b) STM electrochemical cell, c) tweezers, d) glassware, e) glue used for coating of the tips.

The contaminants can be dust, dirt, oil, pigments, rust, grease, algae, fungus, bacteria, lime scale, polishing compounds, flux agents, fingerprints, soot wax and mold release agents, biological soil like blood, etc.

The adopted cleaning procedures involve the use of ultrapure water, bidistilled twice from a permanganate solution (KMnO₄) or piranha solution, a mixture of sulfuric acid (H₂SO₄) and hydrogen peroxide (H₂O₂).

To do this, University of Padua uses:

SONICA®
3200 EP

TECHNICAL CHARACTERISTICS:

Tank dimension: 300X240X100 mm	US Peak Power: 360 W
Capacity Litres: 6	Power Heating: 175 W
External dimension: 325X270X260 mm	Drain cock: Ø 7 mm
Frequency: 40 Khz	Degas Sweep System

Biochemistry laboratories at the “Center for Nanoscience and Technology” of the Italian Institute of Technology in Milano (Italy)

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Anna Paola Lucarelli, PhD and junior Post Doc at the Italian Institute of Technology in Milan - Center of Nanoscience and Technology - describes the study and analysis process performed in the laboratory and explains how the Sonica ultrasonic cleaner is used to extract proteins.

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In the Biochemistry Lab at the CNST-IIT in Milano, researchers focus on the X-ray structure determination of proteins and other biological molecules for a diverse set of applications such as drug design or tissue engineering and for the elucidation of fundamental biological processes.

Natural, correctly folded proteins are produced in large quantities by transforming non-pathogenic bacteria with the gene encoding for the protein of interest. This procedure has been made possible by the development and the optimization of recombinant DNA techniques that had, over the last few decades, an enormous impact on biological research.

Small-scale screenings for optimal bacterial growth conditions are often crucial to improve protein quality and production yield. In fact, many experimental parameters can play a key role in the production process. These include: nature of the medium, its oxygen concentration, nature of the bacterial strain, concentration of the transcription-triggering chemical agent, and others. Such preliminary studies allow for the precise identification of the optimal experimental conditions at which large-scale protein production will need to be carried out.


Figure 1 - Stainless steel holder for 20 vials

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Figure 2 - SDS GEL

Figura 2 - SDS GEL

Following induction and culture growth, bacterial cells are lysed and the protein production pattern is analyzed by gel electrophoresis.

SONICA Ultrasonic bath is an extremely useful tool for all the researchers in the Biochemistry laboratories at the Center for Nanoscience and Technology of the Italian Institute of Technology. The ultrasonic cavitation generated by the Sweep System Technology allows for the simultaneous lysis of multiple bacterial samples of roughly 200 μ -liters volumes each. Furthermore, with the SONICA Ultrasonic bath (digital version) different sonication times and temperature can be conveniently set and maintained constant throughout the entire sonication process, thus allowing for the precise identification of the ideal cell lysis conditions.

Finally, large scale protein production affords sizeable amounts of the protein of interest. At this stage, good quality sonication procedures are key to obtaining the largest possible amounts of correctly folded protein. Further purification steps will then be carried out in order to separate the protein of interest from the endogenous bacterial proteins.

Subsequent functional and structural studies are performed on the purified protein of interest to characterize its intrinsic properties and behavior.

Figure 3 - Sonica 3200 EP S3 ultrasonic bath in the IIT - Italian Institute of Tecnology in Milan.

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“ Such preliminary studies allow for the precise identification of the optimal experimental conditions at which large-scale protein production will need to be carried out.

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“ **The ultrasonic cavitation generated by the Sweep System Technology allows for the simultaneous lysis of multiple...**

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To do this, Institute of Technology in Milan uses:

SONICA®
3200 EP

+ Stainless steel holder for 20 vials

TECHNICAL CHARACTERISTICS:

Tank dimension: 300x240x100 mm	US Peak Power: 360 W
Capacity Litres: 6	Power Heating: 175 W
External dimension: 325x270x260 mm	Drain cock: Ø 7mm
Frequency: 40 Khz	Degas Sweep System

Ultrasonic cleaner for cleaning neutron detector mechanics

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Among the research applications of Sonica cleaners, one of the most interesting is at the Forschungs- Neutronenquelle Heinz Maier-Leibnitz (FRM II). Indeed, the main components of the German research reactor need to be cleaned perfectly. Furthermore, the cleaner is also used to degas the materials that could contaminate the Helium-3 gas.

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The Forschungs-Neutronenquelle Heinz Maier-Leibnitz (FRM II) is the leading German research reactor. It is situated in Garching near Munich and it is operated by the Technische Universität München.

The FRM II serves as neutron source used mainly for fundamental research in physics, chemistry, biology and materials science. More than 25 scientific instruments installed at the FRM II, operated by several universities and institutes, are used to investigate the atomic structures and dynamics of soft and solid matter in a non destructive manner.

Using neutron detectors scientists are able to record scattering patterns of neutrons passing through samples and obtain important insights into the structure and properties of matter.

Main task of the detector & electronics group at the FRM II is developing and maintaining neutron detectors for scientific instruments.

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“ **The 90 liters model can easily accommodate mechanical parts in a size of 400 x 400 mm² allowing a very fast and effective cleaning.**

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
Figure 2 - MWPC Detector: entrance window (Aluminium)

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At present the group is manufacturing a Helium-3 gas-filled Multi-Wire Proportional Counter (MWPC), a two dimensional detector that will be installed at the STRESS-SPEC diffractometer, an instrument devoted to measuring residual stress and texture of materials.

This type of detector provides excellent properties in term of detection efficiency, position resolution and rejection of gamma rays, however in order to achieve and keep these performances some key issues deserve special attention.

An extremely clean atmosphere inside the detector must be guaranteed avoiding the use of any outgassing material that might contaminate Helium-3 gas as well as residuals and dirt have to be accurately removed before closing the detector.

The detector & electronics group has used the Ultrasound Cleaner SONICA 90L EP to clean the MWPC detector vessel and other parts.

The 90 liters model can easily accommodate mechanical parts in a size of 400 x 400 mm² such as stainless steel flange, aluminum window and PCBs allowing a very fast and effective cleaning.


Figure 3 - MWPC Detector: back flange (Steel)

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SONICA®
90L EP

TECHNICAL CHARACTERISTICS:

Tank dimension: 600x500x300 mm	US Peak Power: 2000 W
Capacity Litres: 90	Power Heating: 2000 W
External dimension: 660x560x510 mm	Drain cock: 1/2"
Frequency: 40 Khz	Degas Sweep System

Figure 4 - MWPC Detector: PCB

Figura 4 - Rivelatore: PCB

SONICA 4300ETH S3. Speed and capacity for the Gonzaga Institute in Milan.

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From the most specialized companies to the world of education. The larger Sonica ultrasonic cleaners have proved paramount for washing glassware at school science and chemistry laboratories. At the Gonzaga Institute in Milan, an 18-litre cleaner enables students to wash the instruments used in lessons and experiments both thoroughly and fast.

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Figure 1 - Building Gonzaga Institute

Figura 1 - Edificio Istituto Gonzaga

Figure 2 - Ultrasonic cleaning of an instrument to study the communicating vessels

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Istituto Gonzaga (see Fig.1-3) was founded in 1906 by the religious congregation of the "Brothers of the Christian Schools", founded in 1680 by St. John Baptist de La Salle universal educators' patron. The High Schools of Istituto Gonzaga guarantee a high quality and complete educational path (arts, science, humanities) aimed to that breadth of thought which is a peculiar heritage and legacy of our civilization and of the high school tradition by itself.

At the High Schools' scientific laboratories, we tested the ultrasonic bath SONICA 4300ETH S3 (18 liters).

The main application of this machine is the washing of both laboratory glassware and support tools used for teaching laboratories purposes of both Sciences (Chemistry, Biology and Earth Sciences) and Physics. In particular some experiences of Chemistry, carried out by more classes in a short sequence, have the problem of both a large amount of tools to be washed and different kinds of chemicals to be removed. Moreover, the tools have to be quickly reused by the next classes.

As an example, an experience to study the solubility of different solid compounds with different liquid solvents involves the use of about 90 test tubes for each class using both organic and inorganic solvents and

both hydrophilic and hydrophobic solids compounds. In this case, by using the ultrasonic washing machine we were able to wash more than 60 test tubes in 5 minutes, directly into their test tube rack.

In the case of test tubes, another successful application was the complete removal of both stains and spots of iron hydroxide by coupling the use of ultrasonic washing machine and an acid detergent.

In the laboratory of Physics the complete removal of vegetable oils residues poured into graduated cylinders (experiments to study the motion of bodies) was not an easy task.

By using the ultrasonic washing machine also this operation was successfully accomplished coupling ultrasounds and high temperature (40°C). Finally an interesting application was the washing of an instrument to study the communicating vessels (see Fig. 2). This tool (c.a. 60 years old) is particularly delicate and by hand washing it is virtually impossible to remove the oil residues in the narrow brassy joints. Even in this case by using the ultrasonic washing machine we obtained a satisfactory cleaning of glass vessels and we also cleaned the brassy sections.

Therefore, the final judgment is positive and the main advantages are : a faster washing of the laboratory glassware; the ability to clean not hand washable tools.

Prof. Domenico Fraccalvieri

Science teacher and Laboratory Teaching Assistant in both the laboratories of Chemistry and Physics Istituto Gonzaga, Milano

Figure 3 - Main entrance of Gonzaga Institute (Milan)

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SONICA
4300 ETH

TECHNICAL CHARACTERISTICS:

Tank dimension: 330X300X200 mm	US Peak Power: 600 W
Capacity Litres: 18	Power Heating: 500 W
External dimension: 440X340X410 mm	Drain cock: 1/2"
Frequency: 40 Khz	Sweep System

Ultrasonic precision cleaning of yacht winch components

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SOLTEC provided the hardware for the process of winch cleaning, and offered its technical background and expertise in ultrasonics with an 'on site' technical team at the Desafío Español base. The SONICA Ultrasonic Cleaning system has been able to obtain a precision cleaning on more than 2000 mechanical and hydraulic components restoring them to their original efficiency.

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Harken is synonymous throughout the yachting world with top-quality sailboat deck hardware.

This renowned company, founded in 1967 on the basis of a brilliant invention by Chairman Peter Harken (the first Harken ball bearing block), can now provide all sailboat deck hardware, including winches, needed for America's Cup races.

Cleaning of most of the mechanical components used in Harken winches is done by the SONICA 45L EP S3 ultrasonic cleaner, with Sweep System technology produced by SOLTEC Srl of Milan.

During previous editions of the America's Cup, first in 2003 with **One World Challenge** and thereafter in 2007 with the **Desafío Español 2007 team**, SOLTEC acquired a wide range of in-depth technical skills and developed perfect ultrasonic precision cleaning products, for more than 2000 different mechanical and hydraulic parts.

Figure 1/3/4 - Details with obvious signs of dirt and oxidation before ultrasonic treatment

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Figure 2 - Winch parts before ultrasonic treatment

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The parts, made from technologically advanced materials, such as Delrin, Ertalite, Torton Peak, stainless steels such as 316L, 17-4ph, 13/8 and 2205, aluminum alloys 5083, 6061 and 7000, titanium, carbon fibre and composite epoxy laminates were treated with ultrasound, with incredibly outstanding results.

But let's look in detail at how the SONICA ultrasonic cleaning system works.

During an America's Cup race, the winches have to be in perfect working order, otherwise a mixture of nautical grease, salt, oxidation residues and wear particles can even cause the whole system to seize up.

To prevent this, the whole winch is completely dismantled into its component parts, and each piece is placed in an ultrasonic cleaning bath in a special basket and undergoes an ultrasonic cleaning process that lasts for an average of 30 minutes. The ultrasonic cavitation process generated inside a special alkaline washing liquid specifically designed for this **SONICA ML** application produces micro-cavitation bubbles, which have enough energy to detach all dirt particles from parts being cleaned.

After ultrasonic cleaning, the pieces treated are then reassembled, thus making the part as good as new, and ready for the next race.

Thanks to the experience and advice of Massimiliano Zafaro, Harken Italy Assembly Supervisor, the Italian branch of the Group where Harken winches are designed and produced for the whole world, we have been able to test and re-test the effectiveness of our SONICA ultrasonic cleaning system, which has been successfully applied, not only in recent editions of the America's Cup, but also in countless repairs of mechanical winch components which come to servicing centre to be restored to their original factory condition.

Stefano Guidi
Managing Director
SOLTEC Srl

Figure 5 - Detail with obvious signs of dirt and oxidation before ultrasonic treatment

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Figure 6 - Detail after ultrasonic treatment

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Figures 7/8 - Detail before and after ultrasonic treatment

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
10

11

Figure 9 - Various details after SONICA ultrasonic cleaning with SONICA ML concentrated alkaline detergent

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Figure 10/11 - Detail after SONICA ultrasonic cleaning. Even in spring mechanisms there are no grease or dirt residues. The part is ready to be reassembled after maintenance.

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Stefano Guidi
Managing Director
SOLTEC Srl

Figure 12 - Detail of a crown after SONICA ultrasonic cleaning with SONICA ML concentrated alkaline detergent

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Figure 13 - Ultrasonic cleaning made each part as good as new.

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Figure 14 - Even complex or tiny parts can have all their residues removed.

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SONICA®
45L EP S3

TECHNICAL CHARACTERISTICS:

Tank dimension: 500X300X300 mm	US Peak Power: 1200 W
Capacity Litres: 45	Power Heating: 1000 W
External dimension: 600X340X525 mm	Drain cock: 1/2"
Frequency: 40 Khz	Degas Sweep System

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* Available also in "i" integrated version.

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Technical specifications subject to change.

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Automatic Multifunction System operated by ultrasonics with Thermal Disinfection at a 93°C cycle for the treatment of surgical instruments.

SONICA[®] S.A.M.3L


SONICA® ULTRASONIC BATHS

CLEANING - DEGASSING - DISPERSING
EMULSIFYING - DISSOLVING
HOMOGENISING - PROTEIN EXTRACTION

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Copyright© 1994-2015 Tutti i diritti riservati.
SOLTEC S.r.l