## SOLTEC and IRCCS Mario Negri Institute for Pharmacological Research in Milan A partnership for innovation and participatory research

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At the Mario Negri Institute in Milan, researchers are using increasingly modern equipment and technologies.

Here at SOLTEC, we decided to provide them with our sonication equipment free of charge to help them carry out their research.

Diego Baderna, PhD student and researcher at Laboratory of Environmental Chemistry and Toxicology, uses SONICA Ultrasonic baths to extract pollutants from environmental matrices, while the special SONICA Ultrasonic Extractor is used in two European projects to remove PM10-type particulate matter from membranes.

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In the laboratory of Environmental Chemistry and Toxicology (Department of Environment Health Sciences) of "IRCCS Istituto di Ricerche Farmacologiche Mario Negri" in Milan, researchers are studying the impact of environmental pollutants on human health and on the environmental receptors. The SONICA 3200M ultrasonic bath, currently supplied to the laboratory, is a very useful tool in the preparation of soil leachates according to the ANPA (2000) and APAT (2002 and 2004) guidelines which will be then analyzed by in vitro toxicological assays on human cells or ecotoxicological assays.

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**Figure 2** - Preparation of soil leachates according to the ANPA and APAT.

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The large tank of the instrument allows the simultaneous housing of more reaction chambers, generally consisting of 100 or 250 ml glass flasks, depending on the amount of sample to be extracted. In this way, in addition to the simultaneous processing of multiple samples, it is possible to obtain a comparable extraction efficiency among the samples processed in the same series, ensuring the possibility of comparison. The large tank also allows the insertion of blocks of synthetic ice to cool the temperature of the water in the bath. Using the same technique, the SONICA 3200M was also used to obtain aqueous leachates of earth and rock from mechanized tunneling, comminuted solid waste and ashes while, coupled with organic solvents (ethanol, methanol, acetone: hexane mixture), it was used for extraction of organic pollutants from waste shredding, earth and rock from excavation, especially when the use of high temperature extraction (soxhlet) is not suitable to preserve the integrity of the sample or of the analyzed pollutant (for example PAHs are sensible to temperature > 40°C).

Milan, June 30<sup>th</sup> 2014

**Diego Baderna**

*PhD student and junior researcher  
at Laboratory of Environmental  
Chemistry and Toxicology*

4

3

**Figure 3** - Preparation of soil leachates according to the ANPA and APAT.

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**Figure 4** - Glass flask during ultrasonic extraction with Sonica Bath.

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**The SONICA 3200M ultrasonic bath, currently supplied to the laboratory, is a very useful tool in the preparation of soil leachates**

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**Figure 5** - Glass flask after ultrasonic extraction with Sonic a Bath

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To do this, Istituto Mario Negri uses:

**SONICA® 3200 M S3 and Sonica ultrasonic extractor**

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